

Tango XL

Comfort and stability, also for clients with increased care needs

5100 5705

Tango XL shower toilet chair

The Tango XL shower toilet chair is specially designed for clients who require extra space and support. With its reinforced frame and wider seat, the chair offers comfort and safety for users with a higher body weight—without compromising on manoeuvrability or ease of use.

Featuring foldable arm- and footrests, a hygienic seat, and smooth-rolling castors with brakes, the Tango XL combines practicality with comfort. A reliable solution for care facilities that strive for inclusive care without compromise.

Highly Stable Design. The robust frame, low center of gravity and braked castors ensure maximum stability during transfers and care moments.

Ergonomic Use. Equipped with smart handles and four braked wheels to simplify manoeuvring, ensuring high stability on various surfaces and minimizing the strain on the caregiver's joints.

Reliable & Long-Lasting. Built for daily care and backed by a 60-month warranty, offering a dependable performance.

Antibacterial & Corrosion-Resistant. The antibacterial coating with silver-ion supports hygiene and infection control, while the corrosion-resistant frame with smooth surfaces ensures durability and easy cleaning.

 Discover our tools in action through the Interactive Product videos on our website: www.lopital.com



Technical specifications

Max. load	250 kg / 550 lbs
Product weight	24 kg / 52.9 lbs
Seat width	56 cm / 22 inch
Seat depth	47 cm / 18.5 inch
Seat height	53 cm / 20.9 inch



Suitable for people with very slightly to severely reduced mobility



Optional

6950 1104
Plastic bedpan and
lid.



Options:

- Product no. 6950 1104 Bedpan with lid plastic
- Product no. 6975 8044 Seatbelt 2-point
- Product no. 6975 8045 Seatbelt 4-point
- Product no. 6960 8026 Comfort leg support
- Product no. 7057 0508 Wheel 125 mm Directional

Technical specifications

Total length base frame	100 cm / 39.4 inch
Total width base frame	66 cm / 26.0 inch
Turning radius	115 cm / 45.3 inch



For more information also visit our website: www.lopital.com
Lopital Nederland B.V., Laarakkerweg 9, 5061 JR Oisterwijk
Tel +31 (0)13 52.39.300, E-mail info@lopital.nl, www.lopital.nl

